---

**Tearaway Sheath Introducer**

---

Indications for Use	2
Indikationer for anvendelse	3
Gebruiksindicaties	5
Indications d'emploi	7
Anwendungsgebiete	9
Ενδείξεις χρήσης	11
Felhasználási javallatok	13
Indicazioni per l'uso	15
Bruksindikasjoner	17
Indicações de utilização	19
Indicaciones de uso	21
Indikationer för användning	23
Kullanım Endikasyonları	25

## Tearaway Sheath Introducer

### Indications for Use:

These introducers are used for the percutaneous introduction of diagnostic or therapeutic devices, such as catheters and pacing leads, into the vasculature.

### Intended Patient Population and Group:

For use on both adult and adolescent groups. The specific target population is dependent on the delivered device.

### Device Description:

The tearaway sheath introducer assembly consists of an outer peelable sheath and a dilator. The tear-away sheath has a winged hub to initiate the tear in the sheath.

Kit configuration includes:

**Guidewire:** There are two configurations; Mandrel or Fully Coiled. Mandrel guidewire construction uses a tapered core wire which is soldered or welded to a coil at the distal end. Fully Coiled guidewire construction uses a tightly wound coil which surrounds a tapered core and a safety ribbon.

**Needle:** Introducer needles are composed of two components: a stainless steel cannula with an over molded hub. The introducer needle provides an access path into the vasculature.

### Contraindications:

Use of the introducer is contraindicated if the patient has a known or suspected obstruction in the vessel. There is increased risk of pneumothorax for the patient who has severe chronic lung disease. Poor healing may result in the patient who has had irradiation to the anterior chest.

### Potential Complications:

The potential complications related to the use of the introducer include, but are not limited to the following: Air embolism, device embolism, device dislodgement, pneumothorax, vein thrombosis hematoma formation, hemothorax, vessel erosion trauma to vessels, sepsis.

### Precautions:

Store in a dry, dark, cool place. Do not use if package is open or damaged. Inspect all components prior to use.

### USA CAUTION

**Federal (USA) law restricts this device to sale by or on the order of a physician.**

### Cautions:

- This procedure should only be performed by physicians thoroughly trained in this procedure.
- If resistance is met when advancing or withdrawing the guidewire or the introducer, determine the cause by fluoroscopy and correct before continuing with the procedure.
- Because of the delicate and fragile nature of guidewires, extra care in handling must be taken.
- Guidewires should be routinely inspected prior to use and discarded should any deformities be present in the guidewire.
- Do not use forceps to break the handle and/or to peel the sheath as this may damage the sheath and cause premature withdrawal of the sheath from the patient.
- Do not attempt to use a guidewire over the maximum diameter specified on the package label.
- Individual patient anatomy and physician technique may require procedural variations.
- Insertion into artery may cause excessive bleeding and/or other complications.
- When locking dilator into the outer sheath, Do Not over-tighten dilator. Damage to the locking feature in the sheath will occur if dilator is rotated more than 1/8 turn in the sheath.
- Symmetrical peeling of the sheath is critical.
- Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State.

### Warnings:

- Do not alter this device in any way.
- Do not use alcohol, acetone or solutions containing these agents. These solutions may affect the properties of the plastic components resulting in degradation of the device.
- Do not reuse this device. Reuse will result in increased biocontamination risk for the patient resulting in infection or pyrogenic response.
- Do not attempt to straighten a wire that has been kinked or bent.
- Do not advance a guidewire that is kinked or becomes kinked or bent.
- Do not rotate the guidewire if significant resistance is felt.
- Do not withdraw guidewire through metal needles; guidewire may shear or unravel.
- Do not resterilize.

### USE STERILE TECHNIQUE. A suggested procedure:

1. Peel open package and place contents on sterile field.
2. Prep skin and drape in area of anticipated veni-puncture as desired.
3. Distend the vessel of choice following standard hospital practice for venipuncture. The vein will be much easier to locate if the patient is well hydrated.
4. Insert needle into vessel. The needle position should be verified by observing venous blood return.
5. The angle of the needle should be adjusted depending on the patient's build: shallow in a thin person, deeper in a heavyset person.
6. Aspirate the puncture needle using the syringe.
7. Remove the syringe and insert soft tip of the guidewire through the introducer needle into the vessel. Advance guidewire to required depth. Leave an appropriate amount of guidewire exposed. At no time should the guidewire be advanced or withdrawn when resistance is met. Determine the cause of resistance before proceeding. Fluoroscopic verification of the guidewire position is suggested.
8. Hold guidewire in place and remove introducer needle. Do not withdraw the guidewire back into the cannula as this may result in separation of the guidewire. The cannula should be removed first.
9. Thread the dilator/sheath assembly over the guidewire.
10. Advance the dilator and sheath together with a twisting motion over the guidewire and into the vessel. Fluoroscopic observation may be advisable. Attaching a clamp or hemostat to the proximal end of the guidewire will prevent inadvertently advancing the guidewire entirely into the patient.
11. Remove the vessel dilator and guidewire, leaving the sheath as a conduit into the vessel. Immediately place a finger over the remaining sheath orifice to prevent excessive bleeding or possible air aspiration.
12. Advance the catheter or device through the sheath and into the vessel.  
Note: If desired, lubricate the device with a light coating of appropriate sterile lubricant to facilitate passage through the introducer sheath. Apply moderate pressure until the device tip pushes through the distal end of the sheath. Some resistance may be felt when advancing a device having tines or a flange tip.
13. When the catheter or other device is properly located, peel the sheath away from the device by using the ends of the sheath as handles and pulling them so that the sheath separates longitudinally and breaks at the distal end.

**Caution: Do not attempt to withdraw a tined or flanged device through the introducer sheath. If the device must be withdrawn, remove the sheath first and then gently pull the device out of the vessel. Repeat the insertion procedure with a new introducer.**

**Caution:** When the guidewire is in a vessel, do not advance the movable core if the tip is in a curved shape. Never twist or force the core because excessive force may cause it to penetrate the coil and damage the vessel.

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
**EMERGO EUROPE**  
 Prinsessegracht 20  
 2514 AP The Hague  
 The Netherlands

**Tearaway Sheath Introducer**

<b>Rx ONLY</b>				
Caution: Federal law restricts this device to sale by or on the order of a (licensed healthcare practitioner)	Do Not Use if Package Is Damaged	Single Use Only	Not Made With Natural Rubber Latex	Do Not Re-Sterilize
Ethylene Oxide Sterilized	Instructions for Use	Caution	Non-Pyrogenic	Store In A Dark Place
Keep Dry	Part Number	Lot Number	Manufacturing Date	Use By Date
Manufacturer	Authorized Representative		Quantity	Medical Device
<b>NiTi</b>	<b>NiTi [Au]</b>	<b>NiTi [Pt]</b>		
Nitinol	Nitinol Gold	Nitinol Platinum	Single Sterile Barrier system	Single Sterile Barrier System with Protective Packaging Outside

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**Galt Medical Corp**
U.S.A.