

1

Theragenics

UMBILICAL TAPE

DESCRIPTION

Umbilical Tape is composed of cotton obtained from North American sources and is available in white.

INDICATIONS

Umbilical Tape is indicated for use in tying off the umbilical cord of newborn infants.

ACTIONS

None known.

CONTRAINDICATIONS

None known.

WARNINGS

Do not re-sterilize. Do not re-process. Unit package is sterile unless it has been opened or damaged. Do not use if product packaging is wet. Discard open, unused tape.

In case device is re-used, it is impossible to guarantee its mechanical properties for intended use and there is a risk of cross-infection.

Store at room temperature. Avoid prolonged exposure to elevated temperatures.

Risk of biological contamination; please follow regional or national regulation for elimination of potential contaminated devices, for used tape disposal.

PRECAUTIONS

In handling this product, care should be taken to avoid damage. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

ADVERSE REACTIONS

Do not use this product in patients with known or suspected allergies to this material.

3

2

HOW SUPPLIED

Umbilical tape single use cotton strips are available in sterile, white strands. The tape is supplied sterile in 1/8 inch wide by 18 inches, 30 inches, and 36 inches length in box quantities of 12 and 40 envelopes.

STERILITY

Umbilical Tape is sterilized by Ethylene Oxide.

CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician or licensed practitioner.

DISPOSAL

Risk of biological contamination; please follow regional or national regulation for elimination of potential contaminated devices, for used tape disposal.

4

SYMBOL DEFINITIONS



Lot Number



Expiration Date



Do Not Reuse



Do Not Resterilize



See Instructions For Use



Sterilized By Ethylene Oxide



Keep away from sunlight and heat



Do not use if package is damaged



Keep dry



Manufacturer



CP Medical Inc.
1775 Corporate Drive, Suite 150
Norcross, GA 30093 USA

Theragenics Corporation®

5203 Bristol Industrial Way
Buford, GA 30518 USA
877-444-7333
www.theragenics.com

PUB504 R01

Effective Date 08/23/2022

L